

New! ABC Chair

Provides custom seating for children with special needs

- Offers many points of adjustable support
- Used in school, home and/or rehab setting
- Unique base adapts to unevenness of floors in motion or at rest
- Tilttable seat
- Multiple options allow for custom seating



Trunk supports in various sizes are height and width adjustable and have quick-release fastening


The full-featured, user-friendly ABC chair allows therapists and caregivers to seat each special needs child in an appropriate position, enabling the child to better participate in rehabilitation, home or school activities. The ABC Chair can be customized using the

standard features and optional equipment specific to each child's needs. When properly positioned, a child may experience a decrease in undesirable trunk and leg extensor tone, improved pressure distribution and pelvic and trunk stability.



Call toll free:
1-800-323-5547
In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
In Canada fax 1-905-566-9207



ABC Chair Specifications

Item #	Description
4938	ABC Chair 12"W x 14"D (30.5 x 35.6 cm) complete seat with arm rests and 11"W x 13"H (27.9 x 33 cm) high back. Fits children from ages 3-6.
4939	ABC Chair 14"W x 16"D (35.6 x 40.6 cm) complete seat with arm rests and 12.5"W x 14"H (31.8 x 35.6 cm) high back. Fits children from ages 7-10.
4939A	Hip Abductor
4939B	Lapbelt
4938C	Seat Cover, Incontinent 12 x 14" (30.5 x 35.6 cm)
4939C	Seat Cover, Incontinent 14 x 16" (35.6 x 40.6 cm)
4939KL	Knee Support Left
4939KR	Knee Support Right
4939LL	Legrest, Left
4939LR	Legrest, Right
4939F	One Piece Footboard
4939N	Headrest, Angle Adjustable
4939NSL	Trunk Support, Narrow, Short, Left
4939NSR	Trunk Support, Narrow, Short, Right
4939WLL	Trunk Support Wide, Long, Left
4939WLR	Trunk Support Wide, Long, Right
4939S	Side Support, Long, Each
4939T	Tray

Each ABC Chair includes:

- FLEXBASE underframe with two rear casters which adapt to unevenness in floors when in motion or at rest
- A foot brake to provide stability during transfers or when a locked position is desired
- A manual crank which allows forward or backward tilt
- Gas-spring assisted adjustment to raise or lower seat
- Armrests adjustable for width and height
- Standard adjustable arm rests
- High backrest

Optional equipment includes:

- Trunk supports are adjustable for height, width and depth and have quick-release fastening
- Leg rests which can be individually adjusted for angle, height and depth
- Leg rest with one-piece footboard
- Knee supports which are adjustable for height and depth
- Tray for performance of activities
- Hip abductor
- Lapbelt
- Angle-adjustable headrest



Call toll free:
1-800-323-5547
 In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
 In Canada fax 1-905-566-9207

Available from your local Sammons Preston Rolyan Distributor or call Sammons Preston Rolyan