

Cricket Gait Walker

Adjustable gait trainer for children with motor and balance challenges



Encourages Mobility and Interaction

Indoors and out, the Cricket Walker makes gait training easy and fun for children with motor and balance challenges. The Cricket Walker stimulates walking movements, improves opportunities for social interaction and serves as a link to more advanced walking aids.


Adjustable Walker for Children Ages 2 - 8 Years Old


As attractive as it is functional, the Cricket Walker provides outstanding freedom of movement and exceptional comfort for children ages 2 - 8 years old. Simple height adjustments accommodate a child's growth, and optional accessories provide extra support and comfort for those who need it.



Call toll free:
1-800-323-5547
In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
In Canada fax 1-905-566-9207





The Support Ring (included) swings up and away for easy entry or can swing down to encircle the child for support.

Cricket Walker Specifications

Cricket Walker	Small 1745S	Large 1745L
Maximum Length	31.2" (79.2 cm)	33.9" (86.1 cm)
Maximum Height	28.1" (71.4 cm)	32" (81.3 cm)
Frame Height	10.1" - 14.0" (25.7 - 35.6 cm)	14.0" - 17.9" (35.6 - 45.5 cm)
Frame Height w/Saddle	13.5" - 20.4" (35 - 52 cm)	17.7" - 24" (45 - 62 cm)
Max. User's Height	43" (110 cm)	53" (135 cm)
Max. Weight Capacity	77 lbs. (35 kg)	77 lbs. (35 kg)
Support Ring	8.5" diameter (21.6 cm)	8.5" diameter (21.6 cm)
Ring Height Adjustment	9" - 12.9" (22.9 - 32.8 cm)	9" - 12.9" (22.9 - 32.8 cm)
Weight	16.5 lbs. (7.5 kg)	18.7 lbs. (8.5 kg)

Note: Measure inseam of child with shoes on

Customize the Cricket Walker to Suit a Child's Needs and Level of Independence

- **Curved Saddle** provides extra support while seated
- **Small Trunk Support** keeps smaller children from slouching forward
- **Seat Cushion** provides cushioning between the legs for those who are more mobile, and require little rest time. Can also be placed over saddle for extra seating comfort.
- **Back Support** properly positions child, and provides fall protection



Available Accessories

- 1746S** Saddle
- 1746C** Curved Saddle
- 1746T** Trunk Support
- 1746TS** Small Trunk Support
- 1747-01** Back Support for Small Cricket Walker
- 1747-02** Back Support for Large Cricket Walker
- 1746SC** Seat Cushion
- 1747-03** Grip Rail
- 1747-04** Feet Divider
- 1747-06** Replacement Adjustment Key
- 1747-07** Extra Long Pole (11") for Support Ring
- 1747-08** Pole Support Without Saddle
- 1747-09** Push Handle for A1745S of A1745L



Call toll free:
1-800-323-5547
 In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
 In Canada fax 1-905-566-9207

Available from your local Sammons Preston Rolyan Distributor or call Sammons Preston Rolyan

