

Carrie® Potty Seat by Tumble Forms®... All Grown Up!

Accommodates adolescents to small adults



Maximum user height is 68" and maximum user weight is 130 lbs.

Frame width is 22³/₄", overall height is 53¹/₂", frame length is 33¹/₂"; with footrest 40¹/₂"

Adjustable tray & footrest

Lockable casters

Little kids soon become big kids. That's why our popular Carrie® Potty Seat is available in three sizes; from adolescents to young adults. The above pictured small adult-size seat is specially designed to meet the toileting needs of users with heights ranging 58" to 68" (147-173 cm) and weighing up to 130 pounds (68 kg).

All seats have the same design features that make the Carrie Potty Seat the popular choice in positioning

equipment for toileting needs. The larger seat for adolescents and small adults provides secure positioning while encouraging toileting through increased flexion. The seat is shaped to position legs in slight flexion and abduction to encourage bladder and bowel functions. It can be positioned over 18" institutional toilets, many residential toilets or used as a free-standing commode—so it's ideal for clinical, school or home use.



Call toll free:
1-800-323-5547
In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
In Canada fax 1-905-566-9207

Other Features Include:

- Two comfortable recline positions of 12° or 25°
- Lateral supports and tray help provide trunk stability
- Flexible elastomer splash guard/seat ring protects seat and is removable for easy cleaning
- Built-in anti-thrust seat helps maintain pelvic positioning
- Tumble Forms brand high-density foam has a washable, seamless coating that is urine and stain-resistant, nontoxic and latex free
- Four lockable casters on base facilitates safety with transfers

Carrie® Potty Seat Includes:

- Adjustable tray
- Height adjustable footrest
- Removable splash guard
- Adjustable Headrest
- Hip and chest harness
- 4 position removable footrest
- Removable bucket for use as a commode

2771P Elementary Carrie Potty Seat*

2771J Junior Carrie Potty Seat*

2771SA Small Adult Carrie Potty Seat*

Carrie® Potty Seat Specifications

Item #	Size	Standing Height	Weight of Child
2771P*	Elementary	48" (125 cm)	60 lbs. (27 kg)
2771J*	Junior	58" (147 cm)	100 lbs. (45 kg)
2771SA*	Small Adult	68" (173 cm)	130 lbs. (68 kg)

*Oversized and/or overweight shipping surcharge applies. Call for shipping quote.

Note: may not fit over all toilets in reclined position.



Tray removes for ease in transferring



Splash guard removes for easy cleaning



Includes removable bucket



Call toll free:
1-800-323-5547
 In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
 In Canada fax 1-905-566-9207

*Available from your local Sammons Preston Rolyan
 Distributor or call Sammons Preston Rolyan*