

# Fluidotherapy®

## Dry Heat Therapy for Upper and Lower Extremities

Describes a dry heat therapy that is used to increase circulation and ROM. Also appropriate for pain control and edema reduction. Fluidotherapy units heat and circulate Cellex® particles resulting in a warm massage for treated extremities. All units are electronically programmable. Adjustable controls allow the regulation of temperature (maximum temperature 125°), airflow (particle speed) and pulse mode (pressure oscillations). Each unit has an automatic timer, a wake-up pre-heater mode, and locking casters. UL, CSA approved. Choose from three models:

### *FT110*

A standard single extremity unit for treating the hand, wrist, elbow, ankle or foot. Uses 30 lbs. of Cellex® medium.

**9236-50**

### *K210P*

A large single extremity unit for treating the knee, hand, elbow, foot and ankle. Uses 75 lbs. of Cellex® medium.

**9236-52**

### *Ultra115*

A double extremity unit for treating two hands or two feet independently or simultaneously. Unit is removable from stand for easier lower extremity access. Uses 40 lbs. of Cellex® medium.

**9236-51**

### *Cellex® Medium*

10 lb. bag.

**9236-53**



9236-51



9236-50



9236-52

Fluidotherapy is a registered trademark of the Chattanooga Group, Inc.