

Tumble Forms 2[®] Feeder Seat[®] System Order Form



Please mark clearly with an X the items you wish to purchase. Fax form to Sammons Preston at 800-547-4333 or call in your order to our Pediatric Helpline at 800-499-0285.

Account No. _____ Ship to: _____
 Date: _____ Name: _____
 PO#: _____ Address: _____
 Buyer: _____ City: _____
 Attention: _____ Phone: () _____

Feeder Seat Positioner Features

- Accommodates infants through adolescents
- Shoulder harness slots allow 4" of vertical adjustment of the shoulder straps to accommodate children of different heights
- Contoured interior has a 90° seat-to-back relationship to provide posturally correct seating
- A 45° hip strap and quick-release H-belt help maintain child's position
- Seamless coverings are washable, odor, urine and stain-resistant and nontoxic

Mobile Floor Sitter (includes both the Feeder Seat and the Mobile Base)

_____ Medium, for children up to 48"	4444M	\$728.70
_____ Large for children up to 60"	4542M	\$882.50
_____ Extra-Large for adolescents up to 72"	4542P	\$1,302.50

Mobile Base (may be ordered separately)

_____ Medium Mobile Base	4444MC	\$431.20
_____ Large Mobile Base	4542N	\$426.25
_____ Extra-Large Mobile Base	4542S	\$666.25

Tray for Feeder Seat

_____ Tray for Small/Medium/Large Feeder Seats	2795T	\$881.20
_____ Tray for Extra Large Feeder Seats	2795TXL	\$881.20

Choose One Color

Blue	Red	Teal	Purple	Product	Price
4542BB	4542BR	4542BT	4542BP	Deluxe Floor Sitter Small	\$427.45
4542AB	4542AR	4542AT	4542AP	Deluxe Floor Sitter Medium	\$503.70
4542DB	4542DR	4542DT	4542DP	Deluxe Floor Sitter Large	\$647.45
4542EB	4542ER	4542ET	4542EP	Deluxe Floor Sitter Extra Large	\$1,249.95

Or Mix & Match

Blue	Red	Teal	Purple	Product	Price
2795NB	2795NR	2795NT	2795NP	Small Feeder Seat Positioner	\$306.20
2795MB	2795MR	2795MT	2795MP	Medium Feeder Seat Positioner	\$399.95
2795XB	2795XR	2795XT	2795XP	Large Feeder Seat Positioner	\$518.70
4542WB	4542WR	4542WT	4542WP	Floor Sitter Wedge, fits Small, Medium & Large Feeder Seats	\$274.95
2795YB	2795YR	2795YT	2795YP	Extra Large Feeder Seat Positioner	\$843.70
4542YB	4542YR	4542YT	4542YP	Floor Sitter Wedge, fits Extra Large Feeder Seats	\$596.20

Feeder Seat System Sizing Chart

Dimension	Small Feeder Seat	Medium Feeder Seat	Large Feeder Seat	Extra-Large Feeder Seat
Overall Height	21" (53cm)	27" (68.5cm)	34" (86cm)	41" (104cm)
Outside Width	10½" (26.5cm)	14" (35.5cm)	17" (43cm)	21" (53cm)
Inside Depth	8" (20cm)	9½" (24cm)	13¼" (33cm)	16½" (42cm)
Inner Width	7¾" (18.7cm)	9" (22.8cm)	12½" (31.7cm)	15" (38cm)
Inner Height	18" (46cm)	23" (58.5cm)	29" (73cm)	35" (89cm)
Child's Height	up to 36" (91cm)	up to 48" (122cm)	up to 60" (150cm)	up to 72" (183cm)

Tumble Forms 2[®]

With Anti-Microbial Protection

Odor Resistance, Increased Durability



Deluxe Floor Sitter includes Feeder Seat and Wedge.

