

# The Discovery and Explorer Gait Trainers Order Form

PATTERSON  
MEDICAL

SAMMONS  
PRESTON

Please mark clearly with an X the items you wish to purchase. Fax form to Sammons Preston at 800-547-4333 or call in your order to our Pediatric Helpline 800-499-0285.

Account No. \_\_\_\_\_ Ship to: \_\_\_\_\_  
 Date: \_\_\_\_\_ Name: \_\_\_\_\_  
 PO#: \_\_\_\_\_ Address: \_\_\_\_\_  
 Buyer: \_\_\_\_\_ City: \_\_\_\_\_  
 Attention: \_\_\_\_\_ Phone: (    ) \_\_\_\_\_

## Gait Trainer Features

- Stable wheelbases are easily interchangeable
- Highly adjustable
- Wheelbases feature directional locking and braking casters, bumper guards, and an anti-scissoring foot divider

## Discovery and Explorer Gait Trainers

____ Explorer Size 1	<b>9244-50</b>	<b>\$2,024.95</b>
____ Explorer Size 2	<b>9244-51</b>	<b>\$2,087.45</b>
____ Explorer Size 3	<b>9244-52</b>	<b>\$2,149.95</b>
____ Discovery Size 1	<b>9244-53</b>	<b>\$1,612.45</b>
____ Discovery Size 2	<b>9244-54</b>	<b>\$1,693.70</b>
____ Discovery Size 3	<b>9244-55</b>	<b>\$1,749.95</b>

## Interchangeable Wheel Kits

____ Explorer Wheel Kit (Fits All Sizes)	<b>9244-59</b>	<b>\$898.70</b>
____ Discovery Wheel Kit Size 1	<b>9244-60</b>	<b>\$458.60</b>
____ Discovery Wheel Kit Size 2 & 3	<b>9244-61</b>	<b>\$481.10</b>

## Positioning Accessories

____ Hip Laterals Size 1 & 2	<b>9244-62</b>	<b>\$128.70</b>
____ Hip Laterals Size 3	<b>9244-63</b>	<b>\$134.85</b>
____ Abdominal Pad Size 1 & 2	<b>9244-64</b>	<b>\$146.10</b>
____ Abdominal Pad Size 3	<b>9244-65</b>	<b>\$153.70</b>
____ Arm Support Assembly Size 1 & 2	<b>9244-66</b>	<b>\$368.70</b>
____ Arm Support Assembly Size 3	<b>9244-67</b>	<b>\$384.95</b>
____ Ankle Prompts Size 1	<b>9244-68</b>	<b>\$177.45</b>
____ Ankle Prompts Size 2	<b>9244-69</b>	<b>\$187.45</b>
____ Ankle Prompts Size 3	<b>9244-70</b>	<b>\$196.20</b>
____ Push Handle	<b>9244-71</b>	<b>\$97.45</b>



The Discovery Gait Trainer



The Explorer Gait Trainer

Frame & Accessories Measurements, inches (mm)			
User Age Range	1.5-4 yrs	4-8 yrs	8-12 yrs
User Standing Ht.	33"-41" (850-1050)	41"-50" (1050-1280)	50"-58.5" (1280-1490)
Max User Wt.	66 lbs. (30kg)	88 lbs. (40kg)	110 lbs. (50kg)
Abductor Height (from ground)	11.8"-16.7" (300-425)	12.8"-19.7" (325-500)	18.9"-25.6" (480-650)
Chest Pad Height (from ground)	21"-26" (535-665)	26"-36.5" (665-850)	33.5"-41.7" (850-1060)
Armrest Height	3" (75)	3" (75)	3" (75)
Width Between	6.7"-8.2"	8.2"-11"	11"-12.5"
Hip Laterals	(170-210)	(210-280)	(280-320)
Width Between	5.9"	7"	9"
Chest Lateral	(150)	(180)	(230)
Armrest Angle	+/- 45°	+/- 45°	+/- 45°
Angle-Body Support from Vertical	20°	20°	20°

Discovery – Indoor, inches (mm)			
	Size 1	Size 2	Size 3
Length	35" (900)	39" (990)	47.3" (1200)
Width	30" (770)	30" (770)	30" (770)
Weight	18 lbs. (8kg)	24 lbs. (11kg)	30 lbs. (13.5kg)
Explorer – Outdoor, inches (mm)			
Length	35" (900)	43" (1100)	49" (1250)
Width	26" (660)	26" (660)	26" (660)
Weight	18 lbs. (8kg)	24 lbs. (11kg)	30 lbs. (13.5kg)

Pediatric Helpline (800) 499-0285 Fax (800) 547-4333