



ITO Pushchair Order Form



Please mark clearly with an X the items you wish to purchase. Fax form to Sammons Preston at 800-547-4333 or call in your order to our Pediatric Helpline at 800-499-0285.

Account No. _____ Ship to: _____
 Date: _____ Name: _____
 PO#: _____ Address: _____
 Buyer: _____ City: _____
 Attention: _____ Phone: () _____

Ito Pushchair

Features of the Ito Pushchair:

- Size 1 meets the requirements of ANSI/RESNA WC/Vol 1 section 19 crash testing.
- Light weight (size 1 25 lbs; size 2 27 lbs)
- Adjustable height handle
- Front swivel caster wheels and fixed solid rear wheels

Description	Item#	Price
____ Size 1 Crash Tested Ito	5506-30	\$3,631.30
____ Size 2 Non Crash Tested Ito	5506-50	\$3,631.30
____ Size 1 Headrest	5506-31	\$159.95
____ Size 2 Headrest	5506-51	\$162.45
____ Size 1 Back Extender	0814-35-080	\$248.80
____ Size 2 Back Extender	0814-35-098	\$281.30
____ Size 1 Pelvic Harness	5506-36	\$212.45
____ Size 2 Pelvic Harness	5506-55	\$227.45
____ Lumbar Cushion	5506-37	\$66.30
____ Size 1 Basket	5506-38	\$134.95
____ Size 2 Basket	5506-57	\$134.95
____ Size 1 Sunroof	5506-39	\$392.45
____ Size 2 Sunroof	5506-58	\$396.30
____ Size 1 Rain/wind Cape	5506-40	\$243.80
____ Size 2 Rain/wind Cape	5506-59	\$262.45
____ Size 1 Grip Rail	5506-41	\$147.45
____ Size 2 Grip Rail	5506-60	\$147.45
____ Small Vest Harness	5506-33	\$229.95
____ Medium Vest Harness	5506-34	\$242.45
____ Large Vest Harness	5506-53	\$268.80
____ H-Harness	5506-35	\$127.45
____ X-Harness	5506-54	\$248.80
____ Abductor	5506-45	\$128.80
____ Size 1 Lateral Supports	5506-46	\$207.45
____ Size 2 Lateral Supports	5506-46-01	\$207.45
____ Footstraps	5534-45	\$36.30



	Size 1	Size 2
Seat Width:	8"–11" (20–28cm)	10"–13" (26–34cm)
Seat Depth:	9"–11" (22–28cm)	11"–14" (28–35cm)
Back Height:	20"–24" (50–62cm)	23"–28" (58–72cm)
Height of Foot Plate:	5½"–11" (14–28cm)	8"–14" (20–36cm)
Max Load:	55 lbs. (25kg)	88 lbs. (40kg)
Max Height of User:	47" (120cm)	57" (145cm)
Total Width:	23" (58.5cm)	25" (63.5cm)
Folded Size:	23" x 13" x 30" (58.5 x 33 x 77cm)	25" x 13" x 33" (63.5 x 33 x 83cm)
Weight:	25 lbs. (11kg)	27 lbs. (12kg)