

Leckey® Easy Seat/Potty Trainer/ Waistcoat Order Form



Please mark clearly with an X the items you wish to purchase. Fax form to Sammons Preston at 800-547-4333 or call in your order to our Pediatric Helpline at 800-499-0285.

Account No. _____ Ship to: _____
 Date: _____ Name: _____
 PO#: _____ Address: _____
 Buyer: _____ City: _____
 Attention: _____ Phone: () _____

Features of the Easy Seat

- Attractive, versatile seating solution that grows with the child without compromising positioning needs
- Incorporates seat depth adjustment, leg height adjustment, backrest angle adjustment and height adjustable armrests
- Available in four sizes, use in the home, school, or clinical environments

Easy Seat Complete

Includes tray, 1" side pads, footplate w/straps, 2 swivel/2 fixed casters and lap strap/hip belt.

_____ Size 1	5649-16	\$1,156.20
_____ Size 2	5648-64	\$1,191.20
_____ Size 3	5648-66	\$1,287.45
_____ Size 4	5648-68	\$1,373.70

Easy Seat Basic

Includes non-wheeled legs, hip belt and 1" side pads.

_____ Size 1	5649-17	\$749.95
_____ Size 2	5648-65	\$787.45
_____ Size 3	5648-67	\$874.95
_____ Size 4	5648-69	\$947.45

User Weight & Height

Size	Max. Weight (kg)	Min. Height (mm)	Max Height (mm)
1	121 lbs. (55)	26" (650)	39" (1000)
2	121 lbs. (55)	33" (850)	44" (1120)
3	176 lbs. (80)	39" (1000)	50" (1270)
4	176 lbs. (80)	44" (1120)	59" (1500)

Easy Seat Specifications

Size (mm)	1	2	3	4
Seat Width				
Min.	5" (125)	6" (150)	9" (225)	11" (275)
Max.	9" (225)	10" (250)	13" (325)	15" (375)
Seat Depth				
Min.	7" (175)	9" (225)	10" (250)	12" (300)
Max.	10" (250)	12" (330)	14" (350)	16" (400)
Backrest Height				
Min.	9" (225)	10" (265)	13" (325)	16" (400)
Max.	11" (285)	14" (340)	17" (425)	20" (500)
Seat to Footrest				
Min.	6" (150)	9" (220)	10" (250)	12" (300)
Max.	9" (230)	12" (300)	14" (350)	17" (420)
Seat to Armrest				
Min.	4" (110)	5" (130)	6" (160)	7" (190)
Max.	6" (150)	7" (170)	8" (210)	10" (240)

Tray Dimensions

Size (mm)	1 & 2	3 & 4
Dimensions	18"L x 16.5"W (450 x 420)	23"L x 20.5"W (590 x 520)



Easy Seat Patterns

Meadow Size 1	Mondeo Sizes 2, 3 & 4
-----------------------------	-------------------------------------

The Easy Seat comes with the fabric pattern indicated by the seat size.

Leckey® Easy Seat Order Form (page 2 of 3)



Pommel



Headrest and
headrest laterals supports

Easy Seat Accessories (Size 1)

___ Caster Legs (2 Fixed, 2 Swivel)	9268-75-01	\$246.20
___ Tray	9268-75-02	\$231.20
___ Side Pads ½" (pair)	5649-20	\$77.45
___ Side Pads 1" (pair)	5649-21	\$96.20
___ Side Pads 2" (pair)	5649-22	\$94.95
___ Pommel	5649-23	\$96.20
___ Flexible Lateral Support	5649-24	\$147.45
___ Headrest	5649-25	\$222.45
___ Headrest Lateral Support	5649-26	\$89.95
___ Footplate Raiser (1")	5649-65	\$98.70
___ Footplate w/ Straps	9268-75-11	\$107.45
___ Padded Seat Insert	5649-27	\$72.45
___ Easy Potty Trainer Insert	9243-72	\$63.70

Easy Seat Accessories (Size 2)

___ Caster Legs (2 Fixed, 2 Swivel)	9268-72-01	\$246.20
___ Tray	9268-72-02	\$231.20
___ Side Pads ½" (pair)	5648-75	\$86.20
___ Side Pads 1" (pair)	5648-76	\$104.95
___ Side Pads 2" (pair)	5648-77	\$101.20
___ Pommel	5648-78	\$103.70
___ Flexible Lateral Support	5648-79	\$152.45
___ Headrest	5648-80	\$228.70
___ Headrest Lateral Support	5648-81	\$96.20
___ Footplate w/ Straps	9268-72-11	\$118.70
___ Footplate Raiser (1")	5649-66	\$109.95
___ Padded Seat Insert	5648-72	\$78.70
___ Easy Potty Trainer Insert	9243-75	\$64.95

Easy Seat Accessories (Size 3)

___ Caster Legs (2 Fixed, 2 Swivel)	9268-73-01	\$246.20
___ Tray	9268-73-02	\$231.20
___ Side Pads ½" (pair)	5648-82	\$103.70
___ Side Pads 1" (pair)	5648-83	\$113.70
___ Side Pads 2" (pair)	5648-84	\$119.95
___ Pommel	5648-85	\$114.95
___ Flexible Lateral Support	5648-86	\$159.95
___ Headrest	5648-87	\$239.95
___ Headrest Lateral Support	5648-88	\$107.45
___ Footplate w/ Straps	9268-73-11	\$123.70
___ Footplate Raiser (1")	5649-67	\$118.70
___ Padded Seat Insert	5648-73	\$89.95
___ Easy Potty Trainer Insert	9243-78	\$92.45

Easy Seat Accessories (Size 4)

___ Caster Legs (2 Fixed, 2 Swivel)	9268-74-01	\$246.20
___ Tray	9268-74-02	\$231.20
___ Side Pads ½" (pair)	5648-89	\$113.70
___ Side Pads 1" (pair)	5648-90	\$116.20
___ Side Pads 2" (pair)	5648-91	\$127.45
___ Pommel	5648-85	\$114.95
___ Flexible Lateral Support	5648-86	\$159.95
___ Headrest	5648-87	\$239.95
___ Headrest Lateral Support	5648-92	\$107.45
___ Footplate w/ Straps	9268-74-11	\$131.20
___ Footplate Raiser (1")	5649-68	\$124.95
___ Padded Seat Insert	5648-74	\$98.70
___ Easy Potty Trainer Insert	9243-81	\$92.45

Leckey® Potty Trainer Order Form

Features of the Leckey Potty Trainer

- Molded ABS plastic commode bucket with splashguard, shatter-resistant, crack-resistant, stain-odor resistant
- Backrest height and angle adjustment
- 45° hip belt
- Adjustable side pads
- Height adjustable armrests
- Independently height adjustable legs
- Compatible with all Easy Seat accessories

Leckey Potty Trainer

Includes non-wheeled, adjustable legs, hip belt and 1" side pads

___ Size 1	5649-46	\$749.95
___ Size 2	5649-47	\$781.20
___ Size 3	5649-48	\$912.45
___ Size 4	5649-49	\$962.45



Shown with
optional casters

Leckey® Waistcoat Order Form (Page 3 of 3)

Features of the Waistcoat

- Provides trunk and shoulder positioning
- Two sets of vertical and horizontal straps firmly secure it in position on almost any seating system
- Easily adjusted or quickly removed

Waistcoat.

_____ Size 1	5649-19	\$224.95
_____ Size 2	5648-71	\$224.95
_____ Size 3	9212-53	\$253.70
_____ Size 4	9212-54	\$288.70

Accessories (Size 1)

_____ Replacement Accessory Bib	5649-28	\$59.95
_____ Groin Strap	9212-47	\$56.20

Accessories (Size 2)

_____ Replacement Accessory Bib	5648-70	\$59.95
_____ Groin Strap	9212-48	\$59.95

Accessories (Size 3)

_____ Replacement Accessory Bib	9212-45	\$59.95
_____ Groin Strap	9212-49	\$77.45

Accessories (Size 4)

_____ Replacement Accessory Bib	9212-46	\$59.95
_____ Groin Strap	9212-50	\$93.70



Waistcoat Measurements

Size (mm)	1	2	3	4
Chest Min.	14" (350)	18" (460)	24" (600)	31" (770)
Chest Max.	20" (510)	26" (660)	33" (840)	45" (1100)