

Leckey® Horizon Stander Order Form

Please mark clearly with an X the items you wish to purchase. Fax form to Sammons Preston at 800-547-4333 or call in your order to our Pediatric Helpline at 800-499-0285.

Account No. _____ Ship to: _____
 Date: _____ Name: _____
 PO#: _____ Address: _____
 Buyer: _____ City: _____
 Attention: _____ Phone: () _____

Features of the Leckey Horizon Stander

- Provides appropriate weight bearing for ages 5-18 in prone, supine and vertical position
- Has powered, adjustable actuator to provide effortless movement and positional changes

Horizon Frame. Includes footplate, knee supports, hip plate, chest plate with flexible laterals, battery and charger. (Requires Base and Prone or Supine Pelvic Belt.)

Size 1	Size 2
0815-20-832	0815-20-840
\$2,287.45	\$2,406.20

Horizon Base. Includes four locking casters.

Size 1	Size 2
0815-11-641	0815-11-658
\$2,581.20	\$2,729.95

Accessories

	Size 0	Size 1	Size 2
_____ Prone Pelvic Belt	5505-64 \$134.95	5505-65 \$148.70	5505-66 \$162.45
_____ Prone Tray		5505-62 \$366.20	5505-63 \$424.95
_____ Supine Pelvic Belt	5505-73 \$134.95	5505-74 \$148.70	5505-75 \$162.45
_____ Supine Tray		5505-71 \$298.70	5505-72 \$322.45
_____ Hip Laterals		5505-86 \$243.70	5505-87 \$273.70
_____ Flat Head Support		5505-76 \$194.95	5505-77 \$212.45
_____ Laterals for Flat Head Support		5505-78 \$184.95	5505-78 \$184.95
_____ Contoured Head Support*		5649-03 \$267.45	5649-04 \$267.45
_____ Interface Frame for Contoured Head Support		9208-24-04 \$397.45	9208-24-04 \$397.45
_____ Hyperextension Pads**		5649-69 \$312.45	5649-70 \$338.70
_____ Replacement Battery		5669-99 \$189.95	5669-99 \$189.95
_____ Replacement Battery Charger		5670-98 \$184.95	5670-98 \$184.95



Product Measurement	Size 1	Size 2
L x W x H (horizontal)	44 $\frac{3}{4}$ " L x 22 $\frac{7}{8}$ " W x 30 $\frac{3}{4}$ " H (1135 x 580 x 780 mm)	50 $\frac{1}{2}$ " L x 25 $\frac{1}{4}$ " W x 31 $\frac{3}{4}$ " H (1280 x 640 x 805 mm)
Weight of product	90 lbs. (41kg)	106 lbs. (47kg)
User Standing Height	39 $\frac{1}{2}$ " - 55" (1000-1400 mm)	55" - 71" (1400-1800 mm)
Max User Weight	110 lbs (50kg)	176 lbs (80kg)
Shoulder Height (footplate to top of pad)	32 $\frac{1}{16}$ " - 45 $\frac{1}{16}$ " (820-1150 mm)	39 $\frac{1}{2}$ " - 57" (1000-1450 mm)
Hip Pad Height (footplate to top of pad)	17 $\frac{3}{16}$ " - 29 $\frac{1}{2}$ " (450-750 mm)	29 $\frac{1}{2}$ " - 37 $\frac{1}{2}$ " (750-960 mm)
Knee Height (footplate to middle of pad)	9 $\frac{1}{8}$ " - 17 $\frac{1}{4}$ " (250-450 mm)	15 $\frac{1}{2}$ " - 21 $\frac{1}{2}$ " (400-550 mm)
Width between	8" - 11 $\frac{1}{8}$ "	9 $\frac{1}{2}$ " - 15"
Chest Laterals	(200-300 mm)	(240-380 mm)
Footplate Angle Adjustment (+/-)	10°	10°
Tray Sizes	16 $\frac{1}{2}$ " x 17 $\frac{1}{16}$ " (420 x 440 mm)	20 $\frac{1}{2}$ " x 23 $\frac{1}{4}$ " (520 x 590 mm)
Tray adjustment (Prone)	0° to 40°	0° to 40°
Tray adjustment (Supine)	20°/30°	20°/30°
Caster Size	4" (100mm)	5" (125mm)

*The Contoured Head Support requires an Interface Frame for Head Support (9208-24-04).

**The Hyper Extension Pads require a Contoured Head Support and Interface Frame for Head Support (9208-24-04).

Sandles	Size 0	Size 1	Size 2	Size 3	Size 4
	9244-88	9244-15	9244-25	9244-91	9245-02
	\$132.45	\$134.95	\$141.20	\$153.70	\$153.70