

Leckey® Pronestander Order Form



Please mark clearly with an X the items you wish to purchase. Fax form to Sammons Preston at 800-547-4333 or call in your order to our Pediatric Helpline at 800-499-0285.

Account No. _____ Ship to: _____
 Date: _____ Name: _____
 PO#: _____ Address: _____
 Buyer: _____ City: _____
 Attention: _____ Phone: () _____

The Leckey Pronestander includes activity bowl with tray cover, flexible lateral supports, pelvic band, pelvic plate, knee supports, footplate and sandals.

Pronestander

___ Size 1	5649-37	\$2,218.70
___ Size 2	5648-99	\$2,487.45
___ Size 3	5649-00	\$2,823.70

Accessories (Size 1)

___ Sandals	920821-01-01	\$111.20
___ Pelvic Band	5649-38	\$134.95
___ Sandal Raisers (ea)	9208-21-11	\$54.95
___ Suspension Rail	9208-22-01	\$148.70
___ Knee Spacers (ea)	9208-23-01	\$47.45
___ Contoured Head Support*	5649-40	\$267.45
___ Interface Frame for Head Support	9208-24-04	\$397.45
___ Dynamic Extension Roll	9208-24-05	\$157.45
___ Hyper Extension Pads**	5649-69	\$312.45

Accessories (Size 2)

___ Sandals	920821-01-02	\$123.70
___ Pelvic Band	5648-97	\$148.70
___ Sandal Raisers (ea)	9208-21-12	\$54.95
___ Suspension Rail	9208-22-02	\$148.70
___ Knee Spacers (ea)	9208-23-02	\$47.45
___ Contoured Head Support*	5649-03	\$267.45
___ Interface Frame for Head Support	9208-24-04	\$397.45
___ Dynamic Extension Roll	9208-24-05	\$157.45
___ Hyper Extension Pads**	5649-69	\$312.45

Accessories (Size 3)

___ Sandals	920821-01-03	\$141.20
___ Pelvic Band	5648-98	\$162.45
___ Sandal Raisers (ea)	9208-21-13	\$54.95
___ Suspension Rail	9208-22-03	\$148.70
___ Knee Spacers (ea)	9208-23-03	\$47.45
___ Contoured Head Support*	5649-04	\$267.45
___ Interface Frame for Head Support	9208-24-04	\$397.45
___ Dynamic Extension Roll	9208-24-05	\$157.45
___ Hyper Extension Pads**	5649-70	\$338.70

* The Contoured Head Support also requires item 9208-24-04, Interface Frame for Head Support.

**The Hyper Extension Pads are available only in Mondeo. They require a Contoured Head Support and Interface Frame for Head Support.



Pronestander Patterns

Meadow

Size 1

Mondeo

Sizes 2, 3 & 4

The Pronestander comes with the fabric pattern indicated by the stander size.

Knee Supports

___ Size 1	9208-24-09	\$94.95
___ Size 2	9208-24-10	\$94.95
___ Size 3	9208-24-11	\$94.95
___ Size 4	9208-24-12	\$94.95

Pronestander Specifications

	Size 1	Size 2	Size 3
Chest Plate Height	18.5"-27" (470-680mm)	27"-38" (680-975mm)	38"-53" (975-1350mm)
Base	32" x 24" (820-600mm)	42" x 26" (1070 x 650mm)	50" x 28" (1270 x 700mm)
Activity Tray	16" x 21" (400 x 525mm)	18" x 23" (450 x 575mm)	20" x 25" (500 x 600mm)
User Height	31.5"-41" (800-1050mm)	41"-55" (1050-1400mm)	55"-73" (1400 x 1850mm)
Weight Capacity	132 lbs. (60kg)	176 lbs. (80kg)	198 lbs. (90kg)