



# Leckey Squiggles Saddle Seat Order Form



Please mark clearly with an X the items you wish to purchase. Fax form to Sammons Preston at 800-547-4333 or call in your order to our Pediatric Helpline 800-499-0285.

Account No. \_\_\_\_\_ Ship to: \_\_\_\_\_  
 Date: \_\_\_\_\_ Name: \_\_\_\_\_  
 PO#: \_\_\_\_\_ Address: \_\_\_\_\_  
 Buyer: \_\_\_\_\_ City: \_\_\_\_\_  
 Attention: \_\_\_\_\_ Phone: (     ) \_\_\_\_\_

## Leckey Squiggles Saddle Seat

Benefits include:

- Fully adjustable back rest with lateral supports and chest harness
- 4-point pelvic positioning harness
- Choice of base options offering varying degrees of height adjustment and tilt in space capability
- Lightweight, compact seat can be transferred in seconds from one base to another or removed for transportation and easy storage



**Saddle Seat Shell.** Includes seat base, small 4-point pelvic harness, backrest prompt set and base lock mechanism. (requires seat cover and Hi-Low or Easy Base)

\_\_\_\_ 5624-50 \$1,148.70

**Seat Covers.** Includes: backrest cushion, seat base cushion, upper leg cushion, sacral cushion and pelvic harness.

\_\_\_\_ 5624-51 Green \$308.70

\_\_\_\_ 5624-52 Orange \$308.70

\_\_\_\_ 5624-53 Blue \$308.70

\_\_\_\_ 5624-54 Pink \$308.70

**Hi-Low Base.** Manual gas assist with foot pedal activation.

\_\_\_\_ 5619-75 \$1,511.20

**Easy Base.** With manual crank adjustment.

\_\_\_\_ 5619-76 \$1,371.20

**Push Bar Handle.** Can be used with Hi-Low or Easy Base.

\_\_\_\_ 5619-77 \$218.70

**Footrest Assembly**

\_\_\_\_ 5649-55 \$104.90

**Sandals.** Includes sandal straps. (Requires Footrest Assembly)

\_\_\_\_ 5520-70 Small: 6.9" L x 3"W \$124.95

\_\_\_\_ 5520-71 Medium: 8.3" L x 3.2"W \$134.95

**Sandal Raisers, 1"**

\_\_\_\_ 5535-11 (Single piece fits both sandals) \$68.70

**Rigid Lateral Supports.** Includes hardware and covers.

\_\_\_\_ 5624-62 Green \$256.20

\_\_\_\_ 5624-63 Orange \$256.20

\_\_\_\_ 5624-64 Blue \$256.20

\_\_\_\_ 5624-65 Pink \$256.20

**Chest Support Harness** (Requires rigid lateral supports)

\_\_\_\_ 5619-84 Green \$144.25

\_\_\_\_ 5619-85 Orange \$144.25

\_\_\_\_ 5619-86 Blue \$144.25

\_\_\_\_ 5619-87 Pink \$144.25

**4-point Pelvic Harness, Black, Medium Size**

\_\_\_\_ 5624-55 \$181.20



Pediatric Helpline (800) 499-0285 Fax (800) 547-4333