

Maclaren® Major Elite Push Chair Order Form

PATTERSON
MEDICAL

SAMMONS
PRESTON

Please mark clearly with an X the items you wish to purchase. Fax form to Sammons Preston at 800-547-4333 or call in your order to our Pediatric Helpline at 800-499-0285.

Account No. _____ Ship to: _____
 Date: _____ Name: _____
 PO#: _____ Address: _____
 Buyer: _____ City: _____
 Attention: _____ Phone: () _____

Maclaren Major Elite Push Chair Features

- Swivel-front wheels
- Square-grip handles
- Three-position, detachable foot rests
- Brakes on each rear wheel
- Five-point safety harness
- Hard-wearing, easy-to-clean electric blue upholstery
- Stroller easily folds to 46" x 9¾" x 8¼"
- Accommodates children up to 140 lbs.



Folds like an umbrella stroller



_____ Maclaren Major Elite **Blue** **Charcoal**
 4572 0814-53-042 \$849.95

Optional Accessories

	Blue	Charcoal	
_____ Shopping Basket	0814-53-380	0814-53-315	\$72.45
_____ Padded Seat	0814-53-398	0814-53-067	\$104.95
_____ Sunshade	0814-53-406	0814-53-075	\$97.45
_____ Lateral Inserts (pair)	0814-53-414	0814-53-083	\$62.45
_____ Chest Pad	0814-53-422	0814-53-091	\$47.45
_____ Wind/Rain Shield (must be used with Sunshade)	0814-53-430	0814-53-109	\$81.20
_____ Stroller Boot	0814-53-448	0814-53-117	\$162.45



Shown with
Optional
Accessories

Specifications

Width of seat	15" (38cm)
Depth of seat	10" (25cm)
Height of seat back	20½" (52cm)
Height of seat from ground: front/back	23"/19" (56/48cm)
Overall width	21" (53cm)
Height of handles from the ground	41" (104cm)
Distance from edge seat to foot rest top/middle/bottom	8"/10"/14" (20/26/36cm)
Weight of Stroller: without footrest/with foot rest	14¾ lbs / 16½ lbs 6¾ kg / 7½ kg