



# Soft-Touch™ Therapy Rolls & Wedges Order Form

Please mark clearly with an X the items you wish to purchase. Fax form to Sammons Preston at 800-547-4333 or call in your order to our Pediatric Helpline 800-449-0285.

---

Account No. \_\_\_\_\_ Ship to: \_\_\_\_\_  
 Date: \_\_\_\_\_ Name: \_\_\_\_\_  
 PO#: \_\_\_\_\_ Address: \_\_\_\_\_  
 Buyer: \_\_\_\_\_ City: \_\_\_\_\_  
 Attention: \_\_\_\_\_ Phone: (    ) \_\_\_\_\_

---

### Soft Touch Therapy Rolls feature:

- Integrated Soft-Touch™ Handle for easy control
- Withstands cracking
- Includes anti-microbial agents
- Soft Touch™ Coating

### Soft Touch Therapy Wedges feature:

- Stand On-End™ Design for security and stability
- Soft Touch™ Coating
- Withstands cracking
- Includes anti-microbial agents

#### Soft Touch Rolls.

Chocolate	Diameter/Length	
___ 5533-50	4" x 24"	\$161.20
___ 5533-51	6" x 24"	\$196.20
___ 5533-52	8" x 24"	\$241.20
___ 5533-53	8" x 36"	\$309.95
___ 5533-54	10" x 36"	\$358.70
___ 5533-55	12" x 36"	\$437.45
___ 5533-56	12" x 48"	\$553.70
___ 5533-57	14" x 48"	\$581.20
___ 5533-58	16" x 48"	\$749.95
___ 5535-58	3" x 8"	\$83.70

#### Soft Touch Wedges.

Chocolate	Elevation/Width/Length	Incline	
___ 5533-63	4" x 20" x 22"	10°	\$229.95
___ 5533-64	6" x 20" x 22"	16°	\$256.20
___ 5533-65	8" x 20" x 22"	20°	\$298.70
___ 5533-66	6" x 24" x 28"	12°	\$328.70
___ 5533-67	8" x 24" x 28"	16°	\$358.70
___ 5533-68	10" x 20" x 22"	26°	\$366.20
___ 5533-69	10" x 24" x 28"	20°	\$467.45
___ 5533-70	12" x 24" x 28"	24°	\$478.70



# Soft-Touch™ Sitters

Page 2 of 2

PATTERSON  
MEDICAL

SAMMONS  
PRESTON®

Soft Touch Sitters feature:

- 5-point harness
- Contoured Head Support
- Anti-Thrust Seat Surface
- Up to 25° of tilt with the Tilt Wedge and Mobile Tilt Wedge

five year  
warranty



Base Sizing Chart:	Small		Large	
	Stationary Wood Tilt Wedge	Mobile Wood Tilt Wedge	Stationary Wood Tilt Wedge	Mobile WOOD Tilt Wedge
Use with Sitter Sizes	1,2,3	1,2,3	4,5	4,5
Overall Width - in. (cm)	17.5	17.5	20.5	20.5
Overall Depth - in. (cm)	30	30	35	35
Overall Height - in. (cm)	16	19.5	23	26.5
Weight of Base - lbs. (kg)	12	14	15	17
Degree of Tilt*	25	25	25	25
*Available for Tilt & Mobile Tilt Wedges				

## Special Tomato™ Soft-Touch™ Sitters

Chocolate w/ tan harness	Chocolate w/ lilac harness		
___ 5533-76	5538-76	Size 1 Sitter Only	\$364.95
___ 5533-77	5538-77	Size 2 Sitter Only	\$481.20
___ 5533-78	5538-78	Size 3 Sitter Only	\$611.20
___ 5533-79	5538-79	Size 4 Sitter Only	\$676.20
___ 5533-80	5538-80	Size 5 Sitter Only	\$1029.95
___ 5533-81	5538-81	Size 1 Sitter with Wood Tilt Wedge Kit (No casters)	\$596.20
___ 5533-82	5538-82	Size 2 Sitter with Wood Tilt Wedge Kit (No casters)	\$714.95
___ 5533-83	5538-83	Size 3 Sitter with Wood Tilt Wedge Kit (No casters)	\$862.45
___ 5533-84	5538-84	Size 4 Sitter with Wood Tilt Wedge Kit (No casters)	\$868.70
___ 5533-85	5538-85	Size 5 Sitter with Wood Tilt Wedge Kit (No casters)	\$1,359.95
___ 5533-86	5538-86	Optional Small Wood Tilt Wedge (fits sizes 1-3 sitter)	\$479.95
___ 5533-87	5538-87	Optional Large Wood Tilt Wedge (fits sizes 4-5 sitter)	\$591.20
___ 5533-88	5538-88	Size 1 Sitter Mobile Wood Tilt Wedge Kit	\$811.20
___ 5533-89	5538-89	Size 2 Sitter Mobile Wood Tilt Wedge Kit	\$903.70
___ 5533-90	5538-90	Size 3 Sitter Mobile Wood Tilt Wedge Kit	\$1,048.70
___ 5533-91	5538-91	Size 4 Sitter Mobile Wood Tilt Wedge Kit	\$1,157.45
___ 5533-92	5538-92	Size 5 Sitter Mobile Wood Tilt Wedge Kit	\$1,629.95
___ 5533-93	5538-93	Optional Small Mobile Wood Tilt Wedge (fits sizes 1-3 sitter)	\$587.45
___ 5533-94	5538-94	Optional Large Mobile Wood Tilt Wedge (fits sizes 4-5 sitter)	\$718.70

Sitter Sizing Chart					
Dimensions	Size 1	Size 2	Size 3	Size 4	Size 5
Target Weight - lbs. (kg)	20-40 (9.1-18.1)	30-60 (13.6-27.2)	40-80 (18.1-36.3)	50-110 (22.7-49.9)	90-200 (40.8-90.7)
Target Height - in. (cm)	30-40 (75-100)	35-48 (87.5-120)	40-56 (100-140)	50-65 (125 -162.5)	54-74 (135-185)
Target Age - years	1 to 3	3 to 6	6 to 9	9 to 14	13 to Adult
Inside Seat Depth - in. (cm)	7 (17.5)	9.5 (23.75)	11.5 (28.75)	13.5 (33.75)	16 (40)
Inside Seat Width - in. (cm)	8 (20)	10 (25)	12 (30)	14 (35)	16 (40)
Seat Back Height - in. (cm)	21 (52.5)	26 (65)	29 (72.5)	32.5 (82.5)	39 (97.5)
Floor to Seat Ht. UPRIGHT on <i>Stationary</i> Base	10" (25.4)	10" (25.4)	9.75" (24.8)	11.5" (29.2)	12.5" (31.8)
Floor to Seat Ht. TILTED on <i>Stationary</i> Base	10.5" (26.7)	10.5" (26.7)	10.5 (26.7)	13" (33)	14.25" (36.2)
Floor to Seat Ht. UPRIGHT on <i>Mobile</i> Base	12.5" (31.8)	13.25" (33.7)	13" (33)	14.75" (37.5)	16" (40.6)
Floor to Seat Ht. TILTED on <i>Mobile</i> Base	12.75" (32.4)	13.75 (34.9)	13.75" (34.9)	16" (40.6)	17.75" (45.1)
Overall Seat Depth - in. (cm)	12 (30)	14 (35)	16 (40)	19 (47.5)	21 (52.5)
Overall Seat Width - in. (cm)	12 (30)	14 (35)	16 (40)	18 (45)	20 (50)
Overall Seat Height - in. (cm)	25 (62.5)	29 (72.5)	34 (85)	37.5 (95.25)	44 (110)
Weight of Seat - lbs. (kg)	6 (2.7)	10 (4.5)	13 (5.9)	16 (7.25)	19 (8.6)