

Please mark clearly with an X the items you wish to purchase. Fax form to Sammons Preston at 800-547-4333 or call in your order to our Pediatric Helpline at 800-499-0285.

Account No. \_\_\_\_\_ Ship to: \_\_\_\_\_  
 Date: \_\_\_\_\_ Name: \_\_\_\_\_  
 PO#: \_\_\_\_\_ Address: \_\_\_\_\_  
 Buyer: \_\_\_\_\_ City: \_\_\_\_\_  
 Attention: \_\_\_\_\_ Phone: (    ) \_\_\_\_\_

**Squiggles Stander Support Shell.** Includes: hip and chest support with sternum prompt, chest & hip lateral supports, kneecups, footplate and tray attachment. (Requires covers and pivot base, easel base or mobile base).

\_\_\_\_\_ 5624-30 **\$1,639.95**

**Squiggles Stander Covers.** Includes 2 support harnesses, 2 support cushions, 2 pairs of lateral covers and 1 pair of kneepads.

\_\_\_\_\_ 5624-31 Green **\$462.45**  
 \_\_\_\_\_ 5624-32 Orange **\$462.45**  
 \_\_\_\_\_ 5624-33 Blue **\$462.45**  
 \_\_\_\_\_ 5624-34 Pink **\$462.45**

**Pivot Base.** Gas-operated. Can be operated with child in frame. Great as a lift assist option. Includes four locking casters.

\_\_\_\_\_ 5624-48 **\$1,198.70**

**Easel Base.** Manually adjusted. Cannot be altered with child in frame.

\_\_\_\_\_ 5624-49 **\$616.20**

**Mobile Base.** Manually adjusted, static base. Add Mobile Suspension Kit (5619-34) and 24" or 27" Wheels (5619-35 or 5619-36) for indoor and outdoor mobility.

\_\_\_\_\_ 5619-33 **\$687.45**  
 \_\_\_\_\_ 5619-34 Mobile Suspension Kit **\$513.70**

**Wheels.** Requires the Mobile Suspension Kit (5619-34) and the Mobile Base (5619-33)

\_\_\_\_\_ 5619-35 24" Wheels **\$583.70**  
 \_\_\_\_\_ 5619-36 27" Wheels **\$718.70**

**Grey Tray.**  
 \_\_\_\_\_ 5624-37 **\$309.95**

**Padded Tray Insert for 5624-37.**  
 \_\_\_\_\_ 5535-56 **\$249.95**

**Grab Rail for 5624-37.**  
 \_\_\_\_\_ 5520-79 **\$226.20**

**Transparent Tray.**  
 \_\_\_\_\_ 5619-37 **\$503.70**

**Extra Wide Laterals.**  
 \_\_\_\_\_ 0815-14-702 **\$122.45**

**Headrest Support Shell.** (Requires Flat or Contoured Headrest Cushion)  
 \_\_\_\_\_ 5619-89 **\$172.50**

**Contoured Headrest Cushion.** (Requires Headrest Support Shell 5619-89)

\_\_\_\_\_ 5619-90 Green **\$87.45**  
 \_\_\_\_\_ 5619-91 Orange **\$87.45**  
 \_\_\_\_\_ 5619-92 Blue **\$87.45**  
 \_\_\_\_\_ 5619-93 Pink **\$87.45**

**Flat Headrest Cushion.** (Requires Headrest Support Shell 5619-89)

\_\_\_\_\_ 5619-94 Green **\$87.45**  
 \_\_\_\_\_ 5619-95 Orange **\$87.45**  
 \_\_\_\_\_ 5619-96 Blue **\$87.45**  
 \_\_\_\_\_ 5619-97 Pink **\$87.45**

**Flat Headrest Laterals.** (For use with Flat Headrest Cushion.)


\_\_\_\_\_ 5619-98 Green **\$260.00**  
 \_\_\_\_\_ 5619-98-01 Orange **\$260.00**  
 \_\_\_\_\_ 5619-98-02 Blue **\$260.00**  
 \_\_\_\_\_ 5619-98-03 Pink **\$260.00**

**Square Headrest Bracket.** For headrests with a square post.  
 \_\_\_\_\_ 5535-78 **\$53.70**

**Sandals.** (includes sandal straps)  
 \_\_\_\_\_ 5520-70 Small **\$124.95**  
 \_\_\_\_\_ 5520-71 Medium **\$134.95**

**Sandal Raisers, 1".** (Single piece fits both sandals)  
 \_\_\_\_\_ 5535-11 **\$68.70**

**Replacement Knee Supports**  
 \_\_\_\_\_ 9208-24-09 Size 1 **\$94.95**  
 \_\_\_\_\_ 9208-24-10 Size 2 **\$94.95**  
 \_\_\_\_\_ 9208-24-11 Size 3 **\$94.95**



Sandals are universal to all Squiggles products.



# Leckey Squiggles Stander

## Order Form (Page 2 of 2)



Shown on Mobile Base with Suspension kit, 27" wheels and Transparent Tray.

### Squiggles color swatches



#### Squiggles Stander Specifications:

Age (approx):	1-5
Max. User Weight:	48.4 lbs. (22kg)
User Standing Height:	Min: 29.5" (750mm) Max: 43.5" (1110mm)
Chest Width between laterals:	Min: 6.3" (160mm) Max: 9" (230mm)
Chest Support Height:	Min: 21.7" (550mm) Max: 32" (820mm)
Chest Support:	Angle adjustment +/- 30° Depth adjustment 2" (50mm)
Hip Width between hip guides:	Min: 6.3" (160mm) Max: 9" (230mm)
Hip Support Height from footplate to mid buttocks:	Min: 11.8" (300mm) Max: 20.8" (530mm)
Footplate Angle:	10° Plantarflexion/Dorsiflexion
Sandal Size:	Small 6.9" L x 3"W (175-75mm) Medium 8.3" L x 3.2"W (210-80mm)
Tray Height (Prone):	Min: 21.9" (555mm) Max: 32.5" (825mm)
Tray Height (Supine):	Min: 18.9" (480mm) Max: 32.5" (825mm)
Tray Angle Adjustment:	Prone 30°; Supine 40°

#### Stander Support Frame

Weight:	15.4 lbs. (7kg)
Dimensions:	26.8" L x 15.4" W x 15.4" H (680 x 390 x 390 mm)

#### Pivot Base Dimensions (unfolded)

Weight:	12.1 lbs. (5.5kg)
Dimensions:	32.3" L x 21.7" W x 20.3" H (820 x 550 x 515mm)

#### Pivot Base Dimensions (folded)

Weight:	12.1 lbs. (5.5kg)
Dimensions:	27.6" L x 21.7" W x 11.8" H (700 x 550 x 300mm)

#### Easel Base Dimensions (folded)

Weight:	3.3 lbs. (1.5kg)
Dimensions:	27.2" L x 21.3" W x 2.4" H (690 x 540 x 60mm)

#### Mobile Base (without wheels)

Weight:	13.5 lbs (6.1 kg)
Dimensions:	19.3" W x 29.9" D (490 x 760mm)

#### Mobile Base (with wheels)

Weight:	29.5 lbs (13.4kg)
Dimensions:	29.1" W x 29.5" D (740 x 750mm)