









# Carrie® Seat Sizing Guide

## Seat Dimensions

Carrie Bus Seat		4641B Small	4747B Large	
Carrie Seat	4641S Pre-School	4641M Elementary	4641L Junior	4641SA Small Adult
Inner Seat Height 	22"	26"	30"	33"
Inner Seat Depth 	8"	10"	13"	15 1/2"
Inner Seat Width 	8"	10"	12"	14"
Outside Seat Width 	12"	14 1/4"	16 1/4"	18 1/2"

## Recommended Sizing

	Pre-School	Elementary	Junior	Small Adult
Standing Height 	30" - 38"	38" - 48"	48" - 58"	58" - 68"
Weight 	20 - 40 lbs	30 - 60 lbs	50 - 100 lbs	60 - 130 lbs
Sitting Hip Width 	6" - 8"	8" - 10"	10" - 13"	12" - 15"
Thigh Depth 	8" - 11"	10" - 13"	12" - 15"	14" - 18"